

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\ODED_GONEN\Leah_TMS\Leah_TMS\SCOUT

TA: 8.3 s PAT: Off Voxel size: 2.0x2.0x5.0 mm Rel. SNR: 1.00 SIEMENS: trufi

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	10
Dist. factor	50 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	3.45 ms
TE	1.73 ms
Averages	1
Filter	Distortion Corr.(2D)
Coil elements	A24

Contrast

TD	0 ms
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Table position	H
Table position	0 mm
Inline Composing	Off

System

A24	On
V24	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	L >> R
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Resp. control	Off

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
--------------	----

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Dimension	2D
Reordering	Linear
Asymmetric echo	Off
Bandwidth	673 Hz/Px
Flow comp.	No
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slice-sel.

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\ODED_GONEN\Leah_TMS\Leah_TMS\t1_mpr_sag_1mm_Iso_TxRef_325v

TA: 5:14 PAT: 2 Voxel size: 1.0×1.0×1.0 mm Rel. SNR: 1.00 SIEMENS: tfl

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	176
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2250 ms
TE	3.32 ms
Averages	1
Concatenations	1
Filter	Elliptical filter, B1 filter
Coil elements	A24

Contrast

Magn. preparation	Non-sel. IR
TI	1100 ms
Flip angle	8 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off

Normalize	Off
B1 filter	On
Intensity	Medium
Unfiltered images	On
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Multi-slice mode	Single shot
Series	Ascending
Table position	H
Table position	0 mm
Inline Composing	Off

System

A24	On
V24	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Asymmetric echo	Off
Bandwidth	200 Hz/Px
Flow comp.	No
Echo spacing	7.7 ms
RF pulse type	Fast
Gradient mode	Fast*
Excitation	Non-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\ODED_GONEN\Leah_TMS\Leah_TMS\gre_field_mapping_Scout

TA: 0:40

Voxel size: 1.9×1.9×3.0 mm

Rel. SNR: 1.00

USER: gre_field_mapping_MC

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	15
Dist. factor	200 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	150.0 ms
TE 1	4.00 ms
TE 2	5.00 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	V24

Contrast

MTC	Off
Flip angle	20 deg
Fat suppr.	None
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H

Table position	0 mm
Inline Composing	Off

System

A24	Off
V24	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	240 mm
A >> P	240 mm
F >> H	129 mm

Composing

Sequence

Introduction	On
Dimension	2D
Asymmetric echo	Off
Contrasts	2
Bandwidth	521 Hz/Px
Flow comp.	Yes
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\ODED_GONEN\Leah_TMS\Leah_TMS\gre_field_mapping_5pts

TA: 1:23

Voxel size: 3.8×3.8×3.0 mm

Rel. SNR: 1.00

USER: gre_field_mapping_MC

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	15
Dist. factor	200 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	250 ms
TE 1	4.00 ms
TE 2	5.00 ms
TE 3	6.00 ms
TE 4	8.00 ms
TE 5	12.00 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	V24

Contrast

MTC	Off
Flip angle	20 deg
Fat suppr.	None
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

A24	Off
V24	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	240 mm
A >> P	240 mm
F >> H	129 mm

Composing

Sequence

Introduction	On
Dimension	2D
Asymmetric echo	Off
Contrasts	5
Bandwidth	521 Hz/Px
Flow comp.	Yes
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\ODED_GONEN\Leah_TMS\Leah_TMS\test_sms_mgh_1.6x1.6x1.6mm_MB4-CU

TA: 11:36

PAT: 2

Voxel size: 1.6×1.6×1.6 mm

Rel. SNR: 1.00

USER: ep2d_bold_sms_mgh_v22

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	68
Dist. factor	0 %
Position	L0.3 A37.0 H3.6
Orientation	T > C14.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	1090 ms
TE	26 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	A24

Contrast

MTC	Off
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	600
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	48
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
------------------	-------------

Series

Interleaved

Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

A24	On
V24	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A35.8 F27.8
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	200 mm
! A >> P	200 mm
! F >> H	126 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1628 Hz/Px
Free echo spacing	Off
Echo spacing	0.75 ms
EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
Dummy Scans	6
Dummy Scans	4
SMS Factor	4
RF Clip	0
VERSE Factor	1.00
SMS Shift	2
Kernel Size	3x3
Compression Factor	1.00

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\ODED_GONEN\Leah_TMS\Leah_TMS\test_sms_mgh_1.6x1.6x1.6mm_MB5-CU

TA: 11:58 PAT: 2 Voxel size: 1.6x1.6x1.6 mm Rel. SNR: 1.00 USER: ep2d_bold_sms_mgh_v22

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	75
Dist. factor	0 %
Position	L0.3 A37.0 H3.6
Orientation	T > C14.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	962 ms
TE	26 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	A24

Contrast

MTC	Off
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	700
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	48
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
------------------	-------------

Series

Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

A24	On
V24	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A35.8 F27.8
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	200 mm
! A >> P	200 mm
! F >> H	126 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1628 Hz/Px
Free echo spacing	Off
Echo spacing	0.75 ms
EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
Dummy Scans	6
Dummy Scans	4
SMS Factor	5
RF Clip	0
VERSE Factor	1.00
SMS Shift	2
Kernel Size	3x3
Compression Factor	1.00

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\ODED_GONEN\Leah_TMS\Leah_TMS\PA-test_sms_mgh_1.6x1.6x1.6mm_MB4-CU

TA: 1:36 PAT: 2 Voxel size: 1.6x1.6x1.6 mm Rel. SNR: 1.00 USER: ep2d_bold_sms_mgh_v22

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	68
Dist. factor	0 %
Position	L0.3 A37.0 H3.6
Orientation	T > C14.0
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	1090 ms
TE	26 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	A24

Contrast

MTC	Off
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	50
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	48
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
------------------	-------------

Series

Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

A24	On
V24	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A35.8 F27.8
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	200 mm
! A >> P	200 mm
! F >> H	126 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1628 Hz/Px
Free echo spacing	Off
Echo spacing	0.75 ms
EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
Dummy Scans	6
Dummy Scans	4
SMS Factor	4
RF Clip	0
VERSE Factor	1.00
SMS Shift	2
Kernel Size	3x3
Compression Factor	1.00

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\ODED_GONEN\Leah_TMS\Leah_TMS\PA-test_sms_mgh_1.6x1.6x1.6mm_MB5-CU

TA: 1:32 PAT: 2 Voxel size: 1.6x1.6x1.6 mm Rel. SNR: 1.00 USER: ep2d_bold_sms_mgh_v22

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	75
Dist. factor	0 %
Position	L0.3 A37.0 H3.6
Orientation	T > C14.0
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	962 ms
TE	26 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	A24

Contrast

MTC	Off
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	50
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	48
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Intensity	Weak
Slope	25
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
------------------	-------------

Series

Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

A24	On
V24	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 A35.8 F27.8
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	200 mm
! A >> P	200 mm
! F >> H	126 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1628 Hz/Px
Free echo spacing	Off
Echo spacing	0.75 ms
EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
Dummy Scans	6
Dummy Scans	4
SMS Factor	5
RF Clip	0
VERSE Factor	1.00
SMS Shift	2
Kernel Size	3x3
Compression Factor	1.00

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\ODED_GONEN\Leah_TMS\Leah_TMS\gre_field_mapping_1.6x.16x2

TA: 3:34

Voxel size: 1.6×1.6×2.0 mm

Rel. SNR: 1.00

SIEMENS: gre_field_mapping

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	70
Dist. factor	0 %
Position	L0.3 A37.0 F3.8
Orientation	T > C14.0
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	824.0 ms
TE 1	4.92 ms
TE 2	5.94 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A24

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	None
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H

Table position
Inline Composing

0 mm
Off

System

A24	On
V24	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off

Shim mode Standard

Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto

Adjust volume	
Position	L0.3 A37.0 F3.8
Orientation	T > C14.0
Rotation	90.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	140 mm

Composing

Sequence

Introduction	On
Dimension	2D
Asymmetric echo	Off
Contrasts	2
Bandwidth	337 Hz/Px
Flow comp.	Yes

RF pulse type	Normal
Gradient mode	Normal
RF spoiling	On

Table of contents

\\USER

ODED_GONEN

Leah_TMS

Leah_TMS

SCOUT

t1_mpr_sag_1mm_Iso_TxRef_325v

gre_field_mapping_Scout

gre_field_mapping_5pts

test_sms_mgh_1.6x1.6x1.6mm_MB4-CU

test_sms_mgh_1.6x1.6x1.6mm_MB5-CU

PA-test_sms_mgh_1.6x1.6x1.6mm_MB4-CU

PA-test_sms_mgh_1.6x1.6x1.6mm_MB5-CU

gre_field_mapping_1.6x.16x2